


Final Report

Virucidal Efficacy of Doxmed V-1 on SARS-CoV-2

本測試報告僅供貴公司內部評估，不得作為醫療效果之證明。
The testing report is only for internal evaluation,
not allowed for medical effectiveness advertisement.



2020/12/31

Yi-Ling Lin

Principal Investigator

Infectious Diseases Research and P2/P3 Laboratories,

National Biotechnology Research Park

Nankang, Taipei, Taiwan

Date

本測試報告僅供貴公司內部評估，不得作為醫療效果之證明。
The testing report is only for internal evaluation,
not allowed for medical effectiveness advertisement.

Table of Contents

1. OBJECTIVE	2
2. INTRODUCTION	2
3. MATERIALS and METHODS	2 - 3
4. RESULTS	3 - 4
5. FIGURES and TABLES	5 - 7
6. REFERENCES	8

備註:

1. 本報告係呈現對該送件樣本之測試結果，未經生醫轉譯研究中心書面同意，不得將結果分離使用。
The report presents the test results of the submitted sample. The results should not be used separately without the written consent of the Biomedical Translation Research Center (BioTReC).
2. 該樣品是否適合人類使用不屬於本測試的範圍。
Whether the sample is suitable for human use is outside the scope of this test.
3. 本報告中包含的資訊僅供內部評估和監管審查使用，不得公開用於廣告和商業推廣。
The information contained in this report is only for your internal evaluation and regulatory review, and shall not be made public for advertising and commercial promotion.
4. 未獲中央研究院書面同意前，貴公司不得以任何形式使用中央研究院、中央研究院員工及其所屬單位之名稱、標章、徽章及其他標誌，以使大眾認知中央研究院與貴公司商業發展有任何關連。
Without the written consent of Academia Sinica, your company shall not use any forms of the names, seals, badges, and other marks of Academia Sinica, Academia Sinica employees and their affiliated units to make the public recognize any connection between the Academia Sinica and your company's business development.

1. OBJECTIVE

This project was to evaluate the virucidal efficacy of Doxmed V-1, provided by Houchi Chemical Co., Ltd. (合記化學), on SARS-CoV-2 viruses according to the protocols modified from the peer-reviewed paper [1] and the British standard [2].

2. INTRODUCTION

The antiviral potential was tested by using two cell-based antiviral assays, infection rate determined by immunofluorescent assay and disinfection efficacy determined by median tissue culture infectious dose (TCID₅₀) on African green monkey kidney epithelial Vero-E6 cells. Cytotoxicity of the compounds on these cells were determined by using the CCK8 assay.

3. MATERIALS and METHODS

- 3.1. Vero-E6, ATCC CRL-1586, African green monkey kidney epithelial cells.
- 3.2. SARS-CoV-2, TCDC#4.
- 3.3. Human anti-SARS-CoV-2 N protein antibody (provided by Dr. An-Suei Yang in Genomics Research Center, Academia Sinica).
- 3.4. Cell Counting Kit-8 (CCK8) (Sigma).
- 3.5. Doxmed V-1 (Lot no.: 201907171D) was provided by Houchi Chemical Co., Ltd.
- 3.6. High-content imaging, ImageXpress Micro Imaging XL (Molecular Devices).
- 3.7. The Doxmed V-1, serially diluted with cell culture medium (2% Fetal Bovine Serum + Dulbecco's Modified Eagle Medium), were incubated with SARS-CoV-2 (TCDC#4, 100 PFU) at room temperature for 1 h. The mixtures were then added to Vero-E6 cells for 1 day incubation. Cells were fixed and immunostained with anti-SARS-CoV-2 N protein antibody and goat anti-human IgG-Alexa Fluor 488 (A11013, Invitrogen). Cell nucleus was stained with DAPI (4', 6-Diamidino-2-phenylindole

dihydrochloride, D1306, Invitrogen). To quantify viral infection, images were acquired and analyzed using an ImageXpress Micro XLS Wide field High-Content Analysis System. For cell viability test, Vero-E6 cells were treated with different dilutions of Doxmed V-1 with cell culture medium for 1 day at 37 °C. The cell viability was determined by Cell Counting Kit-8 (CCK-8). 50% inhibition concentration (IC50) and 50% cytotoxic concentration (CC50) were calculated by Prism software. [1]

3.8. The virus titer was determined by virus titration based on median tissue culture infectious dose (TCID₅₀). The Doxmed V-1, serially diluted with cell culture medium, was incubated with SARS-CoV-2 (TCDC#4) at room temperature for 1 hr. The mixtures were then 10-fold serially diluted with cell culture medium and added to Vero E6 cells for 4 day incubation. Cells were fixed with 10% formaldehyde for overnight and stained with 0.5% crystal violet for 20 min. The plates were washed with tap water and scored for infection. The TCID₅₀ was calculated by Reed and Muench Method [3].

3.9. The formula of log reduction value of virus titer is as follows.

Virus reduction for each step is:

$$LRV = \frac{\text{Input virus titer/volume} \times \text{Input volume}}{\text{Output virus titer/volume} \times \text{Output volume}}$$

4. RESULTS

The First Cell-based Antiviral Assay, Infection Rate Determined by Immunofluorescent Assay:

To test the antiviral potential of Doxmed V-1 from Houchi Chemical, SARS-CoV-2 pretreated with serially diluted Doxmed V-1 was added to Vero-E6 cells for incubation. The viral infection was quantified by SARS-CoV-2 N protein expression (**Figure 1A**) and cytotoxicity of the compound was determined by CCK-8 assay. The fluorescent signal was quantified by high-content imaging and calculated the IC50 and CC50 (**Figure 1B**). Doxmed V-1 showed anti-SARS-CoV-2 activity with IC50 at 718.8-fold dilution, however high

cytotoxicity was noted with CC50 at 297.8-fold dilution.

The Second Cell-based Antiviral Assays, Disinfection Efficacy Determined by Median Tissue Culture Infectious Dose (TCID₅₀):

We further tested the disinfectant activity of Doxmed-V1. For untreated virus alone control, the SARS-CoV-2 titer was 10^{8.33} TCID₅₀/ml. For the test group, at 6.25-fold, 12.5-fold, and 25-fold dilution of Doxmed V-1, the virus titer was lower than 10³ TCID₅₀/ml, partly due to the cytotoxicity of Doxmed V-1. At 50-fold dilution of Doxmed V-1, the virus titer was lower than 10² TCID₅₀/ml. At 100-fold, 200-fold, and 400-fold dilution of Doxmed V-1, the virus titers were 10^{3.33}, 10^{4.67} and 10^{5.67} TCID₅₀/ml, respectively (**Figure 2**). The Log-Reduction Value (LRV) of Dexmed-V1 at different dilution is shown in **Table 1**. For cytotoxicity test, different concentrations (100-fold and 300-fold dilution) of Doxmed V-1 were 10-fold serially diluted (10⁰ ~ 10⁻⁷), and added to Vero-E6 cells for 24 hrs. The cytotoxicity was monitored with cytopathic effect (CPE) as shown in **Table 2**.

Overall, at 100, 200 and 400-fold dilution, Doxmed V-1 could inactive SARS-CoV-2 with a LRV at 5, 3.66 and 2.66, respectively. The immunofluorescent assay also showed 200-fold and 400-fold dilution completely blocked SARS-CoV-2 infection. The best inhibition against SARS-CoV-2 was noted at 200~400-fold dilution in this cell-based study condition.

5. FIGURES and TABLES




Figure 1. Anti-SARS-CoV-2 activity of Doxmed V-1 by immunofluorescent assay (IFA). (A) The serially diluted Doxmed V-1 was incubated with 100 PFU SARS-CoV-2 at room temperature for 1 h. The mixtures were then added to Vero-E6 cells for 1 day incubation. Cells were fixed and immunostained with anti-SARS-CoV-2 N protein antibody and goat anti-human IgG-Alexa Fluor 488. (B) SARS-CoV-2 N protein expression levels were acquired and analyzed by using a High-Content Analysis System. Cell viability were determined by CCK-8 assay. 50% inhibition concentration (IC50) and 50% cytotoxic concentration (CC50) were calculated by Prism software.




Figure 2. Disinfection efficacy of Doxmed V-1 was determined by TCID₅₀. The serially diluted Doxmed V-1 was incubated with SARS-CoV-2 (TCDC#4) at room temperature for 1 hr. The mixtures were then 10-fold serially diluted and added to Vero-E6 cells for 4 day incubation. Cells were fixed with 10% formaldehyde and stained with 0.5% crystal violet for 20 min. The plates were washed with tap water and scored for infection. The red boxes indicate wells with cytotoxicity. TCID₅₀ was calculated by Reed and Muench Method [3].

Table 1

	LRV of Doxmed V-1 at different dilution						
Fold of dilution	6.25x	12.5x	25x	50x	100x	200x	400x
Log-reduction value	>5.33	>5.33	>5.33	>6.33	5	3.66	2.66

Table 2

Cytotoxicity of Doxmed V-1									
Fold of dilution		10^0	10^{-1}	10^{-2}	10^{-3}	10^{-4}	10^{-5}	10^{-6}	10^{-7}
V1	100x	+ (0.059%)	+ (71.7%)	- (93.5%)	- (94.5%)	- (93.5%)	- (93.2%)	- (94.2%)	- (94.0%)
	300x	+ (13.29%)	- (90.77%)	- (96.8%)	- (98.1%)	- (96.5%)	- (96.5%)	- (98.8%)	- (99.9%)

The number in brackets indicates the ratio of surviving cells.

6. REFERENCES

1. Liang, J. J., Wei, J. C., Lee, Y. L., Hsu, S. H., Lin, J. J., and Lin, Y. L. 2014. Surfactant-modified nanoclay exhibits an antiviral activity with high potency and broad spectrum. *Jour. of Viro.*, 88: 4218-4228.
2. British standard. Chemical disinfectants and antiseptics. Quantitative suspension test for the evaluation of virucidal activity in the medical area. Test method and requirements (phase 2, step 1). UK. 2019. (BS EN 14476:2013+A2:2019).
3. Reed, L. J., and H. Muench. 1938. A simple method of estimating 50 per cent end-points. *Amer. Jour. Hyg.*, 27: 493-497.

本測試報告僅供 貴公司內部評估，不得作為醫療效果之證明。
The testing report is only for internal evaluation,
not allowed for medical effectiveness advertisement.